



Optos*Advance*TM Adding a New User



Contents



Purpose	1
Add an Account	1
Privilege Descriptions	2
Magnify Glass/Smartzoom Box: Change Settings	2
Modality Preferences.....	3

Purpose

This document will provide guidance on how to add a new user to *OptosAdvance*. Users can be set up as an Admin with full privileges or as a User with certain privileges granted.

You must have administrator privileges to create new users in *OptosAdvance*.

Add an Account

- 1 | Login to *OptosAdvance*
- 2 | Select **Settings**
- 3 | Under **User Management**, select **Accounts**. Existing *OptosAdvance* user accounts are shown
- 4 | Select **Add**
- 5 | In the **Account** area, enter the required information, indicated by a yellow font color
 - **User Name:** Username to login to *OptosAdvance*
 - **Role:** Admin, User, or Guest
 - By default:
 - Guests have no privileges
 - Users have basic privileges, such as accessing the Patient Study Directory
 - Administrators have full privileges. Only Administrators can manage users
 - **Email:** Email address
 - **Last Name, First Name, Middle Name, Prefix, Suffix**
 - **Password:** Password to login to *OptosAdvance*
 - **Expiry Date:** Date the user’s access to *OptosAdvance* will expire
 - Select  and select an expiry date. Select whether the user’s account will be locked or deleted on the expiry date.
 - Select  to remove the expiry date and set the user’s access to **Unlimited**



Optos*Advance*TM Adding a New User



- 6 | Review and update user privileges in the **Privileges** section (see [Privilege descriptions](#)). These privileges are initially based on the user's role and group but can be modified.
 - o **Grant a privilege to a user:** Select a privilege in the **Revoked** area, then select **Grant**.
 - o **Remove a privilege from a user:** Select a privilege in the **Granted** area, then select **Revoke**.

Privilege Descriptions

The following privileges can be assigned to a group or user.

BookmarkSaveSend: enables ability to "Capture" image

Collaboration: enables ability to "Bookmark" a specific layout/view which is then saved in the presentation toolbar

Content download: enables ability to download image/series as a DICOM and "Download Study"

Content Upload: enables ability to "Upload New Study" and "Attach files to this study/patient"

CreateAccounts: enables ability to "Add, Edit, or Delete" user accounts

Dicom Delete: enables ability to delete image or delete series

DICOM Print: Mange DICOM Printers

DICOM Store: enables "send series" to Montage

EditPatientStudy: enables "Edit Image Header" for changing laterality and "Edit Patient/Study" for editing patient information as well as "Merge Studies" or "Merge Patients"


MprProtocols: View Slab views.

ThreeDProtocols: View 3D views.

** Note: A full list of Privilege descriptions can be found in the Optos*Advance* Help Menu by searching Privilege descriptions

Magnify Glass/Smartzoom Box: Change Settings

- 1 | Login as the new user
- 2 | Click on **Settings>Preferences>User Preferences**
- 3 | Click on **Magnify Glass**
- 4 | Change Zoom Factor 1X (It defaults to 4X)





Optos*Advance*TM Adding a New User



Modality Preferences

- 1 | From the Settings icon, select Preferences
- 2 | Select **Modality Preferences**
- 3 | Change Modality to **OP**
 - a. Ensure Series Navigation is **None**
 - b. Ensure Create Combines Series is **Yes**
- 4 | Click Save
- 5 | Change Modality to **OPT**
 - a. Ensure Series Navigation is **None**
 - b. Ensure Create Combined Series is **Yes**
 - c. Change Force Isotropic Display to **Yes**
- 6 | Click Save

Note: Modality Preferences must be changed for all users

