



# Optos *Advance*™ Uploading Non-DICOM Images



## Contents

Purpose .....	1
Uploading a Non-DICOM Image.....	1

## Purpose


This document provides guidance on how to upload a non-DICOM image into *OptosAdvance*.

## Uploading a Non-DICOM Image


- 1 | From within *OptosAdvance*, select **Upload** from the toolbar.



- 2 | Select **Upload non-DICOM files**.



- 3 | Drag the non-DICOM files from your system into the rectangle or select Choose files.



Note: You can upload multiple photos of one patient to *OptosAdvance* by selecting them all for import. For multiple patients, this process must be done individually.



# Optos *Advance*<sup>TM</sup> Uploading Non-DICOM Images



4 | Select **Next** after choosing the non-DICOM file(s).

Optos *ADVANCE*<sup>TM</sup>

Drag and drop Non-DICOM files below

Selected 1 Files

Or

Select Non-DICOM files using the button below

Choose Files Choroidal... a-1326206.tif

Back Next Cancel

5 | Select **Reconcile**.

Optos *ADVANCE*<sup>TM</sup>

Enter Patient Study Details

Last Name	
First Name	
Middle Name	
Date Of Birth	yyyyMMdd
Sex	
MRN	
Accession Number	
Study Id	
Study Description	
Series Number	1001
Series Description	

Reconcile... Back Upload Cancel

